PULPDENT®

Embrace™ WetBond™ Resin Cement

Dual Cure • Fluoride Releasing • Hydrophilic • Radiopaque

Embrace Resin Cement is a dual cure, self-adhesive cement that is moisture friendly and forms chemical bonds to dentin and enamel, precious and non-precious metals, ceramics, composites, and fiber posts. Bonding agents are not required if retention form is adequate; however, they may be used if desired. It is not necessary to etch dentin, but etching uncut enamel surfaces is indicated. The low 12-micron film thickness is ideally suited for luting applications. The cement is radiopaque and is filled with both sub-micron and micron sized particles.

Apply Embrace Resin Cement to slightly moist teeth. Remove surface water, but do not desiccate the tooth, which naturally contains water. Ceramic, metal, resin desensitized, prehybridized and cured resin surfaces should be dry. Place an oxygen inhibiting gel over exposed margins to prevent the formation of an oxygen-inhibited layer, if desired.

PRODUCT DESCRIPTION

Dual cure resin cement

INDICATIONS

Recommended for cementation of PFM, gold, CEREC® and reinforced all-ceramic crowns; inlays; onlays; gold, metal, titanium and fiber posts; and for bonding stainless steel and nylon splinting materials.

CONTRAINDICATIONS

Not recommended for porcelain veneers.

SETTING CHARACTERISTICS

Embrace WetBond Resin Cement is dual cure. The self-cure polymerization occurs by anaerobic polymerization. This means that the reaction does not begin while the material is exposed to air. It begins when the restoration is seated and air is eliminated. This feature provides very long working time and explains why material left on the mixing pad may not polymerize.

DIRECTIONS FOR USE

1. Your laboratory should etch or abrade the tooth contacting surfaces of restorations. If not, etch or microabrade these surfaces. Rinse and dry. Restorations should be closely fitted.
2. Rinse tooth preparation. Lightly dry and remove surface water with compressed air or a cotton pellet. Do not desiccate the tooth, which naturally contains water. Tooth surfaces should be shiny or glossy. Overly wet surfaces will result in decreased bond strength. Ceramic, metal, resin desensitized, prehybridized, and cured resin surfaces should be dry.
3. Etching and bonding agents are not required on dentin when there is adequate retention form. Etching uncut enamel is indicated. Bonding posts is recommended.
4. If cementing to existing ceramic, metal, resin desensitized, prehybridized, and cured composite surfaces, clean and etch or mechanically roughen the surface, rinse and dry. Be sure to dry these surfaces before applying Embrace Resin Cement.
5. Place lubricant on adjacent teeth.
6. Automix Embrace Resin Cement by placing a mixing tip on the double barrel syringe and dispensing material (the first time using the syringe discard the first 2-3 mm, which may not have an equal mix of base and catalyst), or dispense equal amount of base and catalyst and hand mix. Carefully recap syringes. Do not cross-contaminate base and catalyst.
7. Place cement inside the crown and seat crown, or place cement directly onto the prepared tooth and seat inlay, onlay or post
8. Maintain positive pressure on restoration for 2 minutes. During this time remove excess from margins with a brush, or tack cure margins for 1-2 seconds and gently tease away excess with a suitable instrument. Clear contacts and interproximal areas with floss. Cement may bond to adjacent teeth if excess is not removed.
9. Apply oxygen inhibiting gel on exposed margins to eliminate the air-inhibited layer, if desired.
10. Light cure all-ceramic restorations and margins of PFM crowns for 20-30 seconds. Light cure posts for 40 seconds. Complete anaerobic self-cure setting time at mouth temperature is 3 minutes.

Note: Apply disposable barrier sleeves/wraps over multiple-use dental dispensers before use with each patient. For additional information, refer to: <https://www.fda.gov/medical-devices/dental-devices/multiple-use-dental-dispenser-devices>

CAUTION

Uncured material may cause eye or skin irritation on contact. Dental professionals should wear safety glasses and surgical gloves. Do not exceed manufacturer's recommended curing time for the light you are using.

STORAGE AND HANDLING

- Store tightly sealed in original container at cool room temperature. Avoid direct light, extremes of temperature, contamination and sources of ignition.
- Shelf life of unopened product: 2 years from date of manufacture.
- Re-cap immediately after use.

REPORTING OF SERIOUS INCIDENTS

Any serious incidents involving the device should be reported to the manufacturer and the Competent Authority of the Member State in which the incident occurred, or other appropriate regulatory body. Reports of serious incidents should be sent to:

Pulpdent Corporation

80 Oakland Street, Watertown, MA 02472 USA

Tel: 617-926-6666 / Fax: 617-926-6262

email: pulpdent@pulpdent.com

Note: Refer to the Safety Data Sheet for a list of materials present in this product.

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Hofmeister Dental B.V.
Lifco Dental AB
UNIDENT AB
PULPDENT® Corporation
80 Oakland Street • Watertown, MA 02472 • U.S.A.